
## **Ultra High Performance Liquid Chromatography (UHPLC)**

HPLC is used frequently in biochemistry and analytical chemistry to separate, identify, and quantify compounds based on their idiosyncratic polarities and interactions with the column's stationary phase. Thermo Scientific $^{\text{TM}}$  Vanquish $^{\text{TM}}$  UHPLC system offer the separation of compounds within shorter time frame than normal HPLC. It comes with UV/Vis detector and charged aerosol detector are available for different application.



For enquiries, please contact the staff in-charge: Ms Maria Chong

Phone: 6790-4851 | Email: MASCHONG@ntu.edu.sg

**Usage Rate Refer <u>Here</u>**